

# Anti-IL-2 Monoclonal Antibody,Unconjugated, Detector (PBS Only)

## AGMB06907

### Description

---

This Anti-IL-2 Monoclonal Antibody,Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB06907
<b>Contents:</b>	100µg Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	Interleukin-2,T-cell growth factor (TCGF)
<b>Applications:</b>	<b>ELISA</b>
<b>Research Area:</b>	Immunology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Rat
<b>Clone:</b>	-
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P17108
<b>Immunogen:</b>	-
<b>Gene ID:</b>	116562
<b>Gene Name:</b>	IL2

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Rabbit

**Isotype:** -

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 2-8°C. Avoid freezing.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 0.01M PBS, pH 7.2

**Antibody Dilution Ratio:** Optimize antibody dilution ratio according to assay requirements and sample type.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.