

PPP2R1A/PPP2R2A/PPP2R1B Antibody

CAB18571

Description

This PPP2R1A/PPP2R2A/PPP2R1B Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	CAB18571
Contents:	20 μ L, 100 μ L Bradford Reagent: 1 vial (2ml)
Category:	Polyclonal Antibody
Synonyms:	PPP2R1A/PPP2R2A/PPP2R1B
Clone:	-
Applications:	WB ELISA
Conjugation:	-
Reactivity:	Human, Mouse, Rat


Antibody Data

Gene ID:	5518 5520 5519
Uniprot:	AB_2862335
Host Species:	Rabbit
Purification:	Affinity purification
Observed MW:	55kDa/66kDa
Calculated MW:	-

Preparation & Storage

Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3. Store Bradford Reagent at Room Temperature for 1 Year.	
Positive Sample:	HeLa, HepG2, K-562, NIH/3T3, C6, Mouse liver, Mouse lung, Rat brain, Rat testis	
Recommended Dilutions:	WB	1:500 - 1:2000
	ELISA	Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol	

Validation Data



Western blot analysis of various lysates using PPP2R1A/PPP2R2A/PPP2R1B Rabbit pAb (CAB18571) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (AbGn00020). Exposure time: 10s.