

SARS-CoV-2 Spike RBD (N501Y) Monoclonal Antibody

CAB21254

Description

This SARS-CoV-2 Spike RBD (N501Y) Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	CAB21254
Contents:	20 μ L, 100 μ L Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Synonyms:	spike glycoprotein, SARS-CoV-2 Spike RBD (N501Y)
Clone:	ARC52673
Applications:	WB ELISA
Conjugation:	Unconjugated
Reactivity:	SARS-CoV-2

Antibody Data

Gene ID:	43740568
Uniprot:	-
Host Species:	Rabbit
Purification:	Affinity purification
Observed MW:	36kDa (Recombinant Protein)
Calculated MW:	141KD

Preparation & Storage


Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS containing 50% glycerol and 0.05% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.
Store Bradford Reagent at Room Temperature for 1 Year.

Positive Sample: Recombinant SARS-CoV-2 Spike RBD (N501Y) Protein(RP01258)

Recommended Dilutions:	WB	1:10000 - 1:40000
	ELISA	Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data



Western blot analysis of recombinant SARS-CoV-2 Spike RBD (N501Y) Protein, using SARS-CoV-2 Spike RBD (N501Y) Rabbit mAb (CAB21254) at 1:10000 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 100 ng per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (AbGn00020). Exposure time: 3s.