

FceR1 alpha Monoclonal Antibody

CAB22729

Description

This FceR1 alpha Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	CAB22729
Contents:	20 μ L, 100 μ L Bradford Reagent: 1 vial (2ml)
Category:	
Synonyms:	-
Clone:	-
Applications:	[App]
Conjugation:	-
Reactivity:	-

Antibody Data

Gene ID:	-
Uniprot:	-
Host Species:	
Purification:	-
Observed MW:	-
Calculated MW:	-

Preparation & Storage

Storage:	- Store Bradford Reagent at Room Temperature for 1 Year.
Positive Sample:	-
Recommended Dilutions:	N/A
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol

Validation Data

