

## NDUFA5 Antibody

CAB3976

### Description

---

This NDUFA5 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                      |                                                           |
|----------------------|-----------------------------------------------------------|
| <b>SKU:</b>          | CAB3976                                                   |
| <b>Contents:</b>     | 20 $\mu$ L, 100 $\mu$ L<br>Bradford Reagent: 1 vial (2ml) |
| <b>Category:</b>     | Polyclonal Antibody                                       |
| <b>Synonyms:</b>     | B13, NUFM, UQOR13, CI-13kB, CI-13KD-B, NDUFA5             |
| <b>Clone:</b>        | -                                                         |
| <b>Applications:</b> | <span>WB</span> <span>ELISA</span>                        |
| <b>Conjugation:</b>  | Unconjugated                                              |
| <b>Reactivity:</b>   | Human, Mouse, Rat                                         |

### Antibody Data


---

|                       |                       |
|-----------------------|-----------------------|
| <b>Gene ID:</b>       | 4698                  |
| <b>Uniprot:</b>       | AB_2765424            |
| <b>Host Species:</b>  | Rabbit                |
| <b>Purification:</b>  | Affinity purification |
| <b>Observed MW:</b>   | 13kDa                 |
| <b>Calculated MW:</b> | 13kDa                 |

## Preparation & Storage

|                                           |                                                                                                                                                                                                                                                        |                                                                                                                             |
|-------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Storage:</b>                           | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.<br>Store Bradford Reagent at Room Temperature for 1 Year.            |                                                                                                                             |
| <b>Positive Sample:</b>                   | HepG2, 293T, NCI-H460, HeLa, Jurkat, Mouse brain, Mouse kidney, Mouse heart, Rat brain                                                                                                                                                                 |                                                                                                                             |
| <b>Recommended Dilutions:</b>             | <b>WB</b>                                                                                                                                                                                                                                              | 1:500 - 1:2000                                                                                                              |
|                                           | <b>ELISA</b>                                                                                                                                                                                                                                           | Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Protein Quantification (Optional):</b> | To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <a href="https://www.assaygenie.com/bradford-protein-assay-protocol/">https://www.assaygenie.com/bradford-protein-assay-protocol/</a> to view the full protocol |                                                                                                                             |

## Validation Data



Western blot analysis of various lysates using NDUFA5 Rabbit pAb (CAB3976) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (AbGn00020). Exposure time: 60s.