

## Mouse anti TAP-Tag Monoclonal Antibody

CABE021

### Description

---

This Mouse anti TAP-Tag Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	CABE021
<b>Contents:</b>	20 $\mu$ L, 100 $\mu$ L Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	Monoclonal Antibody
<b>Synonyms:</b>	TAP, TAP tag, TAP-tag
<b>Clone:</b>	AMC0510
<b>Applications:</b>	<span>WB</span> <span>ELISA</span>
<b>Conjugation:</b>	Unconjugated
<b>Reactivity:</b>	Species independent

### Antibody Data

---

<b>Gene ID:</b>	-
<b>Uniprot:</b>	AB_2770411
<b>Host Species:</b>	Mouse
<b>Purification:</b>	Affinity purification
<b>Observed MW:</b>	75kDa (Recombinant TAP protein)
<b>Calculated MW:</b>	-

## Preparation & Storage

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.


Store Bradford Reagent at Room Temperature for 1 Year.

**Positive Sample:** Recombinant TAP protein

<b>Recommended Dilutions:</b>	<b>WB</b>	1:5000 - 1:10000
	<b>ELISA</b>	Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data



Western blot analysis of recombinant TAP protein using mouse anti TAP-Tag mAb (CABE021) at dilution of 1:5000. Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (CABS003) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST.