

Affinity Gel-conjugated Mouse anti DDDDK-Tag mAb

CABE061

Description

This Affinity Gel-conjugated Mouse anti DDDDK-Tag mAb is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	CABE061
Contents:	20 μ L, 100 μ L Bradford Reagent: 1 vial (2ml)
Category:	Monoclonal Antibody
Synonyms:	DDDDK, DDDDK tag, DDDDK-tag
Clone:	AMC0382
Applications:	
Conjugation:	Affinity Gel
Reactivity:	Species independent


Antibody Data

Gene ID:	-
Uniprot:	AB_2863789
Host Species:	Mouse
Purification:	Affinity purification
Observed MW:	Refer to Figures
Calculated MW:	-


Preparation & Storage

Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3. Store Bradford Reagent at Room Temperature for 1 Year.
Positive Sample:	-
Recommended Dilutions:	N/A
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol

Validation Data



Immunoprecipitation analysis of 300ug extract cell lysate from 293T cells transfected with GSK3B expression vector containing a DDDDK-tag with 30uL Affinity Gel-conjugated Mouse anti DDDDK-Tag mAb {Anti-DDDDK Affinity Gel} (CABE061). Western blot was performed from the immunoprecipitate using Rabbit anti DDDDK-Tag antibody at 1:10000 dilution.



Immunoprecipitation analysis of 300ug extract cell lysate from 293T cells transfected with expression vector containing a DDDDK-tag with 30uL Affinity Gel-conjugated Mouse anti DDDDK-Tag mAb antibody (CABE061). Western blot was performed from the immunoprecipitate using Rabbit anti DDDDK-Tag antibody at 1:10000 dilution.