

## Mouse anti Strep II-Tag mAb

CABE066

### Description

---

This Mouse anti Strep II-Tag mAb is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                      |                                                 |
|----------------------|-------------------------------------------------|
| <b>SKU:</b>          | CABE066                                         |
| <b>Contents:</b>     | 20 µL, 100 µL<br>Bradford Reagent: 1 vial (2ml) |
| <b>Category:</b>     | Monoclonal Antibody                             |
| <b>Synonyms:</b>     | Strep II, Strep II tag, Strep II-tag            |
| <b>Clone:</b>        | AMC0521                                         |
| <b>Applications:</b> | <span>WB</span> <span>ELISA</span>              |
| <b>Conjugation:</b>  | -                                               |
| <b>Reactivity:</b>   | Species independent                             |

### Antibody Data

---

|                       |                       |
|-----------------------|-----------------------|
| <b>Gene ID:</b>       | -                     |
| <b>Uniprot:</b>       | AB_2863792            |
| <b>Host Species:</b>  | Mouse                 |
| <b>Purification:</b>  | Affinity purification |
| <b>Observed MW:</b>   | 32kDa                 |
| <b>Calculated MW:</b> | -                     |

## Preparation & Storage

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Store Bradford Reagent at Room Temperature for 1 Year.


**Positive Sample:** 293T transfected with Strep II-tagged protein (Mouse)

**Recommended Dilutions:**

|              |                                                                                                                             |
|--------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>WB</b>    | 1:2000 - 1:5000                                                                                                             |
| <b>ELISA</b> | Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data



Western blot analysis of lysates from wild type (WT) and 293T cells transfected with Strep II-tagged protein using Mouse anti Strep II-Tag mAb (CABE066). Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (CABS003) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (AbGn00020). Exposure time: 10 s.