

## Phospho-Camk2-T286 Monoclonal Antibody

CABP1386

### Description

---

This Phospho-Camk2-T286 Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	CABP1386
<b>Contents:</b>	20 $\mu$ L, 100 $\mu$ L Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	Monoclonal Antibody
<b>Synonyms:</b>	CAMKA, MRD53, MRT63, CaMKIIalpha, CaMKIINalpha, CAM2, CAMK2, CAMKB, MRD54, CaMKIIbeta, CAMKD, CAMK, CAMKG, MRD59, CAMK-II, Phospho-Camk2-T286
<b>Clone:</b>	ARC57709
<b>Applications:</b>	<span style="background-color: #e91e63; color: white; padding: 2px;">WB</span> <span style="background-color: #2196f3; color: white; padding: 2px;">ELISA</span>
<b>Conjugation:</b>	Unconjugated
<b>Reactivity:</b>	Mouse

### Antibody Data

---

<b>Gene ID:</b>	816 817 818 815
<b>Uniprot:</b>	-
<b>Host Species:</b>	Rabbit
<b>Purification:</b>	Affinity purification
<b>Observed MW:</b>	50 kDa
<b>Calculated MW:</b>	54 kDa

## Preparation & Storage


**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS containing 50% glycerol and 0.05% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Positive Sample:** Mouse brain

<b>Recommended Dilutions:</b>	<b>WB</b>	1:1000 - 1:2000
	<b>ELISA</b>	Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data



Western blot analysis of lysates from Mouse brain using Phospho-Camk2- Rabbit mAb (CABP1386) at 1:1000 dilution incubated at room temperature for 1.5 hours. Mouse brain treated with CIP (1U/ul) at 37°C for 1 hours. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 30 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (AbGn00020). Exposure time: 45 s.