

HRP-conjugated Donkey anti-Goat IgG (H+L)

CABS031

Description

This HRP-conjugated Donkey anti-Goat IgG (H+L) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	CABS031
Contents:	20 µL, 100 µL Bradford Reagent: 1 vial (2ml)
Category:	
Synonyms:	-
Clone:	-
Applications:	WB IHC-P
Conjugation:	HRP
Reactivity:	Goat

Antibody Data

Gene ID:	-
Uniprot:	AB_2769846
Host Species:	Goat
Purification:	Affinity purification
Observed MW:	-
Calculated MW:	-

Preparation & Storage

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH 7.3.
Store Bradford Reagent at Room Temperature for 1 Year.

Positive Sample: -

Recommended Dilutions:	WB	1:5000 - 1:100000
	IHC-P	1:500 - 1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Full page HTML saved