

Bavituximab Antibody

GXBS00471

Description

This Bavituximab Antibody is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXBS00471
Contents:	Bavituximab Antibody (100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	Phosphatidylserine
Form:	Liquid
Storage Buffer:	0.01M PBS pH 7.4
Applications:	ELISA, FACS, Functional assay, Research in vivo

Antibody Data

Host species:	Chimeric
Clonality:	Monoclonal
Accession:	1446756-47-3
Isotype:	IgG1-kappa
Purification:	Protein A/G purified from cell culture supernatant.
Endotoxin Level:	< 10 EU/mg

Preparation & Storage

Stability and Storage:	Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only. Not suitable for clinical or therapeutic use.