

Blinatumomab Antibody

GXBS00574

Description

This Blinatumomab Antibody is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXBS00574
Contents:	Blinatumomab Antibody (100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	B-lymphocyte antigen CD19
Form:	Liquid
Storage Buffer:	30 mM Histidine, pH5.8, 10% sucrose, 0.02% Tween80
Applications:	ELISA, FCM, Functional assay, Research in vivo

Antibody Data

Reactivity:	Human
Host species:	Mouse
Clonality:	Monoclonal
Accession:	P15391 & P07766
Isotype:	(scFv-kappa-heavy)-(scFv-heavy-kappa)
Purification:	Purified by Nickel column.
Endotoxin Level:	< 10 EU/mg

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only. Not suitable for clinical or therapeutic use.