

Efdamrofusp Alfa (Research-Grade)

GXBS01198

Description

This Efdamrofusp Alfa (Research-Grade) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXBS01198
Contents:	Efdamrofusp Alfa (Research-Grade) (100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	C3
Form:	Liquid
Storage Buffer:	0.01M PBS pH 7.4.
Applications:	ELISA, Functional assay, Research in vivo

Antibody Data

Reactivity:	Human
Clonality:	Monoclonal
Accession:	P01024 & P0C0L5 & P15692

Isotype: Fusion - FLT1 (fms-related tyrosine kinase 1, vascular endothelial growth factor receptor 1, VEGFR-1, VEGFR, FLT, FRT, vascular permeability factor receptor) domain 2 fragment (14-204) - KDR (kinase insert domain receptor, vascular endothelial growth factor receptor 2, VEGFR2, VEGF-R2, FLK1, CD309) domain 3 fragment (206-308) - IGHG1 Fc (Fragment constant (205-431)) - CR1 (complement receptor type 1, C3b/C4b receptor, C3BR, Knops blood group, KN, CD35) fragment (1-192)

Purification: Protein A/G purified from cell culture supernatant.

Endotoxin Level: < 10 EU/mg

Preparation & Storage

Stability and Storage: Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Store Bradford Reagent at Room Temperature for 1 Year.

Shipping: Gel Pack

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Note: For research use only. Not suitable for clinical or therapeutic use.