

Bemiltenase Alfa (Research-Grade)

GXBS02225

Description

This Bemiltenase Alfa (Research-Grade) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXBS02225
Contents:	Bemiltenase Alfa (Research-Grade) (100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	Coagulation factor X-activating enzyme heavy chain
Form:	Liquid
Storage Buffer:	0.01M PBS pH 7.4
Applications:	ELISA, Functional assay, Research in vivo

Antibody Data

Reactivity:	Daboia siamensis (Eastern Russel's viper) (Daboia russelii siamensis)
Accession:	Q7LZ61 & Q4PRD1 & Q4PRD2
Isotype:	Daboia siamensis metalloproteinase (Covalent heterotrimer with alpha heavy chain domain and a cysteine-rich region disulfide-bridged to beta light chain disulfide-bridged to gamma light chain).
Purification:	Purified by Nickel column.
Endotoxin Level:	< 10 EU/mg

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only. Not suitable for clinical or therapeutic use.