

Amtabafusp alfa (Research-Grade)

GXBS02518

Description

This Amtabafusp alfa (Research-Grade) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXBS02518
Contents:	Amtabafusp alfa (Research-Grade) (100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	CD3E&gp120/Glycoprotein 120
Form:	Liquid
Applications:	ELISA, FACS, Functional assay, Research in vivo

Antibody Data

Reactivity:	Human immunodeficiency virus type 1 (HIV-1)
Host species:	Human
Clonality:	Monoclonal
Accession:	P07766&Q03805
Isotype:	half-IG G1 - lambda / CD4 -h-CH2-CH3
Endotoxin Level:	< 10 EU/mg

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only. Not suitable for clinical or therapeutic use.