

Arumakimig Antibody

GXBS02592

Description

This Arumakimig Antibody is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXBS02592
Contents:	Arumakimig Antibody (100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	IL18/IL1F4 & IL1B/IL1F2
Form:	Liquid
Applications:	ELISA, FACS, Functional assay, Research in vivo

Antibody Data

Reactivity:	Human
Host species:	Human
Clonality:	Monoclonal
Accession:	Q14116 & P01584
Isotype:	IgG1, lambda2-IgG1, Kappa
Endotoxin Level:	< 10 EU/mg

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only. Not suitable for clinical or therapeutic use.