

Anafiltamig (Research-Grade)

GXBS02640

Description

This Anafiltamig (Research-Grade) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXBS02640
Contents:	Anafiltamig (Research-Grade) (100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	CD79β/CD79B & CD3E & CD20/MS4A1
Storage Buffer:	0.01M PBS pH 7.4
Applications:	ELISA, FACS, Functional assay

Antibody Data

Reactivity:	Human
Host species:	Human
Accession:	P40259 & P07766 & P11836
Isotype:	immunoglobulin (H-gamma1_L-kappa)_(scFvkh-G1(h-CH2- CH3)-20mer-scFvkh)
Endotoxin Level:	< 10 EU/mg

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only.