

Emaretamig (Research-Grade)

GXBS02651

Description

This Emaretamig (Research-Grade) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXBS02651
Contents:	Emaretamig (Research-Grade) (100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	CLDN18.2 & CD3E
Storage Buffer:	0.01M PBS pH 7.4
Applications:	ELISA, FACS, Functional assay

Antibody Data

Reactivity:	Human
Host species:	Human
Accession:	P56856 & P07766
Isotype:	immunoglobulin (H-gamma1_L-kappa)_(G1(VH-CH1-h)-10mer-scFvIh-10mer-G1(h-CH2-CH3)_L-kappa)_x0002_
Endotoxin Level:	< 10 EU/mg

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only.