

Emzotamig (Research-Grade)

GXBS02652

Description

This Emzotamig (Research-Grade) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXBS02652
Contents:	Emzotamig (Research-Grade) (100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	CD335/NCR1/NKp46 & CD269/TNFRSF17/BCMA
Storage Buffer:	0.01M PBS pH 7.4
Applications:	ELISA, FACS, Functional assay

Antibody Data

Reactivity:	Human
Host species:	Human
Accession:	O76036 & Q02223
Isotype:	immunoglobulin [VH-(H-gamma1)_ V-lambda-10mer-V-kappa-10mer-C_x0002_kappa)]_G1(h-CH2-CH3)
Endotoxin Level:	< 10 EU/mg

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only.