

Glarivibart Antibody

GXBS02660

Description

This Glarivibart Antibody is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXBS02660
Contents:	Glarivibart Antibody (100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	F
Form:	Liquid
Storage Buffer:	0.01M PBS pH 7.4.
Applications:	ELISA, FCM, Functional assay

Antibody Data

Reactivity:	Human respiratory syncytial virus A (HRSV-A), Human respiratory syncytial virus B (HRSV-B)
Host species:	Human
Clonality:	Monoclonal
Accession:	P03420, P11209, 12568, O36634, P13843
Isotype:	IgG1, kappa
Endotoxin Level:	< 10 EU/mg

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only.