

Lumivatamig (Research-Grade)

GXBS02671

Description

This Lumivatamig (Research-Grade) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXBS02671
Contents:	Lumivatamig (Research-Grade) (100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	CLDN18.2 & CD3DE & CD3D
Storage Buffer:	0.01M PBS pH 7.4
Applications:	ELISA, FACS, Functional assay

Antibody Data

Reactivity:	Human
Host species:	Human
Accession:	P56856 & P07766 & P04234
Isotype:	immunoglobulin (H-gamma1_L-kappa)_(H-gamma1-17mer_x0002_scFvLh_L-kappa
Endotoxin Level:	< 10 EU/mg

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only.