

Ordastobart Antibody

GXBS02680

Description

This Ordastobart Antibody is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXBS02680
Contents:	Ordastobart Antibody (100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	ACT35 antigen
Form:	Liquid
Storage Buffer:	0.01M PBS pH 7.4.
Applications:	ELISA, FCM, Functional assay

Antibody Data

Reactivity:	Human
Host species:	Human
Clonality:	Monoclonal
Accession:	P43489
Isotype:	immunoglobulin only-heavy-chain G1 (VH-VH-VH-h-CH2-CH3- CHS) dimer
Endotoxin Level:	< 10 EU/mg

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only.