

## Laronidase (Research-Grade)

**GXBS02751**

### Description

---

This Laronidase (Research-Grade) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                      |                                                                              |
|----------------------|------------------------------------------------------------------------------|
| <b>SKU:</b>          | GXBS02751                                                                    |
| <b>Contents:</b>     | Laronidase (Research-Grade) (100 µg, 1 mg)<br>Bradford Reagent: 1 vial (2ml) |
| <b>Form:</b>         | Lyophilized                                                                  |
| <b>Applications:</b> | Functional assay                                                             |

### Antibody Data

---

|                         |                            |
|-------------------------|----------------------------|
| <b>Purification:</b>    | Purified by Nickel column. |
| <b>Endotoxin Level:</b> | < 10 EU/mg                 |

### Preparation & Storage

---

|                                           |                                                                                                                                                                                                                                                         |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Stability and Storage:</b>             | Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.<br>Store Bradford Reagent at Room Temperature for 1 Year.                                                                                         |
| <b>Shipping:</b>                          | Gel Pack                                                                                                                                                                                                                                                |
| <b>Protein Quantification (Optional):</b> | To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <a href="https://www.assaygenie.com/bradford-protein-assay-protocol/">https://www.assaygenie.com/bradford-protein-assay-protocol/</a> to view the full protocol. |

**Note:** For research use only.