

## Zolbetuximab Neutralizing Antibody

SKU: GXDA00019 | Target: Zolbetuximab | Monoclonal | Anti-Drug Antibody

Zolbetuximab Neutralizing Antibody (SKU: GXDA00019) is a monoclonal anti-drug antibody (ADA) supplied by Assay Genie for research into immunogenicity against Zolbetuximab. It is raised in rabbit and is intended for the detection and characterisation of anti-Zolbetuximab antibody responses in ELISA. This clone is characterised as neutralizing and inhibits the interaction of Zolbetuximab with its target, which suits it to neutralising-antibody assay controls and potency work.

Anti-drug antibodies are a central readout in biotherapeutic development. An immune response raised against a molecule such as Zolbetuximab can alter its clearance, reduce exposure and blunt clinical response, so well-characterised ADA reagents are required as assay standards, positive controls and calibrators throughout preclinical and clinical immunogenicity testing.

### Specifications

|                 |                                    |
|-----------------|------------------------------------|
| Product Name    | Zolbetuximab Neutralizing Antibody |
| SKU             | GXDA00019                          |
| Target Drug     | Zolbetuximab                       |
| Antibody Type   | Monoclonal                         |
| Clone ID        | SAA2837                            |
| Host Species    | Rabbit                             |
| Isotype         | IgG                                |
| Applications    | ELISA                              |
| Available Sizes | 50 µg, 100 µg, 1 mg                |
| CAS Number      | 1496553-00-4                       |

### Properties

|              |                                                     |
|--------------|-----------------------------------------------------|
| Format       | Liquid                                              |
| Buffer       | 0.01M PBS pH 7.4                                    |
| Purity       | >95% as determined by SDS-PAGE.                     |
| Purification | Protein A/G purified from cell culture supernatant. |
| Guarantee    | 12 months from date of dispatch                     |

### Storage & Stability

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

### Notes

For research use only. Not for use in diagnostic or therapeutic procedures.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | [www.assaygenie.com](http://www.assaygenie.com) | [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com)  
Datasheet: <https://www.assaygenie.com/content/Antibody%20Genie/GXDA/GXDA00019.pdf> | Version 1.0, 28 July 2026