

Obinutuzumab Non-Neutralizing Antibody

Description

This Obinutuzumab Non-Neutralizing Antibody is supplied as a kit for advanced applications.

Product Information

SKU:	GXDA00049
Contents:	Obinutuzumab Non-Neutralizing Antibody (liquid)
Clone:	SAA2867
Applications:	ELISA

Antibody Data

CAS Number:	949142-50-1
Size:	50 µg, 100 µg, 1 mg
Host Species:	Rabbit
Clonality:	Monoclonal
Isotype:	IgG
Immunogen:	Obinutuzumab
Conjugation:	Unconjugated
Buffer:	0.01M PBS pH 7.4
Purity:	>95% as determined by SDS-PAGE.
Purification:	Protein A/G purified from cell culture supernatant.
Endotoxin Level:	< 10EU/mg

Preparation & Storage

Storage:	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Shipping:	Gel Pack
Preparation:	Recommended dilution: ELISA: 1:2000-1:20000

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.