

Nivolumab Non-Neutralizing Antibody (SAA2349)

SKU: GXDA00221 | Target: Nivolumab | Monoclonal | Anti-Drug Antibody

Nivolumab Non-Neutralizing Antibody (SAA2349) (SKU: GXDA00221) is a monoclonal anti-drug antibody (ADA) supplied by Assay Genie for research into immunogenicity against Nivolumab. It is raised in rabbit and is intended for the detection and characterisation of anti-Nivolumab antibody responses in ELISA. This clone is characterised as non-neutralizing and binds Nivolumab without blocking its interaction with its target, which suits it to total-ADA and bridging assay formats.

Anti-drug antibodies are a central readout in biotherapeutic development. An immune response raised against a molecule such as Nivolumab can alter its clearance, reduce exposure and blunt clinical response, so well-characterised ADA reagents are required as assay standards, positive controls and calibrators throughout preclinical and clinical immunogenicity testing.

Specifications

Product Name	Nivolumab Non-Neutralizing Antibody (SAA2349)
SKU	GXDA00221
Target Drug	Nivolumab
Antibody Type	Monoclonal
Clone ID	SAA2349
Host Species	Rabbit
Isotype	IgG
Applications	ELISA
Available Sizes	50 µg, 100 µg

Properties

Format	Liquid
Buffer	0.01M PBS pH 7.4
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin Level	Please contact with the lab for this information.
Guarantee	12 months from date of dispatch

Storage & Stability

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Background

ELISA: 1:2000-1:20000

Notes

For research use only. Not for use in diagnostic or therapeutic procedures.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
Datasheet: <https://www.assaygenie.com/content/Antibody%20Genie/GXDA/GXDA00221.pdf> | Version 1.0, 28 July 2026