

Adalimumab Non-Neutralizing Antibody (SAA2490)

SKU: GXDA00225 | Target: Adalimumab | Monoclonal | Anti-Drug Antibody

Adalimumab Non-Neutralizing Antibody (SAA2490) (SKU: GXDA00225) is a monoclonal anti-drug antibody (ADA) supplied by Assay Genie for research into immunogenicity against Adalimumab. It is raised in rabbit and is intended for the detection and characterisation of anti-Adalimumab antibody responses in ELISA. This clone is characterised as non-neutralizing and binds Adalimumab without blocking its interaction with its target, which suits it to total-ADA and bridging assay formats.

Anti-drug antibodies are a central readout in biotherapeutic development. An immune response raised against a molecule such as Adalimumab can alter its clearance, reduce exposure and blunt clinical response, so well-characterised ADA reagents are required as assay standards, positive controls and calibrators throughout preclinical and clinical immunogenicity testing.

Specifications

Product Name	Adalimumab Non-Neutralizing Antibody (SAA2490)
SKU	GXDA00225
Target Drug	Adalimumab
Antibody Type	Monoclonal
Clone ID	SAA2490
Host Species	Rabbit
Isotype	IgG
Applications	ELISA
Available Sizes	50 µg, 100 µg
CAS Number	331731-18-1

Properties

Format	Liquid
Buffer	0.01M PBS pH 7.4
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin Level	Please contact with the lab for this information.
Guarantee	12 months from date of dispatch

Storage & Stability

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Background

ELISA: 1:2000-1:20000

Notes

For research use only. Not for use in diagnostic or therapeutic procedures.