

## Vedolizumab Polyclonal Antibody

SKU: GXDA00291 | Target: Vedolizumab | Polyclonal | Anti-Drug Antibody

Vedolizumab Polyclonal Antibody (SKU: GXDA00291) is a polyclonal anti-drug antibody (ADA) supplied by Assay Genie for research into immunogenicity against Vedolizumab. It is raised in rabbit and is intended for the detection and characterisation of anti-Vedolizumab antibody responses in ELISA, IHC, WB.

Anti-drug antibodies are a central readout in biotherapeutic development. An immune response raised against a molecule such as Vedolizumab can alter its clearance, reduce exposure and blunt clinical response, so well-characterised ADA reagents are required as assay standards, positive controls and calibrators throughout preclinical and clinical immunogenicity testing.

### Specifications

|                 |                                 |
|-----------------|---------------------------------|
| Product Name    | Vedolizumab Polyclonal Antibody |
| SKU             | GXDA00291                       |
| Target Drug     | Vedolizumab                     |
| Antibody Type   | Polyclonal                      |
| Host Species    | Rabbit                          |
| Isotype         | IgG                             |
| Applications    | ELISA, IHC, WB                  |
| Available Sizes | 50 µg, 100 µg                   |
| CAS Number      | 943609-66-3                     |

### Properties

|                 |                                                     |
|-----------------|-----------------------------------------------------|
| Format          | Liquid                                              |
| Buffer          | 0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300. |
| Purification    | Purified by antigen affinity column.                |
| Endotoxin Level | Please contact with the lab for this information.   |
| Guarantee       | 12 months from date of dispatch                     |

### Storage & Stability

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

### Background

ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000

### Notes

For research use only. Not for use in diagnostic or therapeutic procedures.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | [www.assaygenie.com](http://www.assaygenie.com) | [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com)  
Datasheet: <https://www.assaygenie.com/content/Antibody%20Genie/GXDA/GXDA00291.pdf> | Version 1.0, 28 July 2026