

Teriparatide Polyclonal Antibody

SKU: GXDA00292 | Target: Teriparatide | Polyclonal | Anti-Drug Antibody

Teriparatide Polyclonal Antibody (SKU: GXDA00292) is a polyclonal anti-drug antibody (ADA) supplied by Assay Genie for research into immunogenicity against Teriparatide. It is raised in rabbit and is intended for the detection and characterisation of anti-Teriparatide antibody responses in ELISA, IHC, WB.

Anti-drug antibodies are a central readout in biotherapeutic development. An immune response raised against a molecule such as Teriparatide can alter its clearance, reduce exposure and blunt clinical response, so well-characterised ADA reagents are required as assay standards, positive controls and calibrators throughout preclinical and clinical immunogenicity testing.

Specifications

Product Name	Teriparatide Polyclonal Antibody
SKU	GXDA00292
Target Drug	Teriparatide
Antibody Type	Polyclonal
Host Species	Rabbit
Reactivity	Human
Isotype	IgG
Applications	ELISA, IHC, WB
Available Sizes	50 µg, 100 µg
CAS Number	52232-67-4

Properties

Format	Liquid
Buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Endotoxin Level	Please contact with the lab for this information.
Guarantee	12 months from date of dispatch

Storage & Stability

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Background

ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000

Notes

For research use only. Not for use in diagnostic or therapeutic procedures.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
Datasheet: <https://www.assaygenie.com/content/Antibody%20Genie/GXDA/GXDA00292.pdf> | Version 1.0, 28 July 2026