

ADI-1942 Antibody (SAA2807)

SKU: GXDA00295 | Target: ADI-1942 | Monoclonal | Anti-Drug Antibody

ADI-1942 Antibody (SAA2807) (SKU: GXDA00295) is a monoclonal anti-drug antibody (ADA) supplied by Assay Genie for research into immunogenicity against ADI-1942. It is intended for the detection and characterisation of anti-ADI-1942 antibody responses in ELISA.

Anti-drug antibodies are a central readout in biotherapeutic development. An immune response raised against a molecule such as ADI-1942 can alter its clearance, reduce exposure and blunt clinical response, so well-characterised ADA reagents are required as assay standards, positive controls and calibrators throughout preclinical and clinical immunogenicity testing.

Specifications

Product Name	ADI-1942 Antibody (SAA2807)
SKU	GXDA00295
Target Drug	ADI-1942
Antibody Type	Monoclonal
Isotype	IgG
Applications	ELISA
Available Sizes	50 µg, 100 µg
Accession	ADI-1942

Properties

Format	Liquid
Buffer	0.01M PBS pH 7.4
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin Level	Please contact with the lab for this information.
Guarantee	12 months from date of dispatch

Storage & Stability

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Background

ELISA: 1:2000-1:20000

Notes

For research use only. Not for use in diagnostic or therapeutic procedures.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
Datasheet: <https://www.assaygenie.com/content/Antibody%20Genie/GXDA/GXDA00295.pdf> | Version 1.0, 28 July 2026