

Anti-REGN5381 Antibody (REGN9035)

SKU: GXDA00302 | Target: REGN5381 | Monoclonal | Anti-Drug Antibody

Anti-REGN5381 Antibody (REGN9035) (SKU: GXDA00302) is a monoclonal anti-drug antibody (ADA) supplied by Assay Genie for research into immunogenicity against REGN5381. It is raised in human and is intended for the detection and characterisation of anti-REGN5381 antibody responses in ELISA, Bioactivity: FACS, Functional assay, Research in vivo.

Anti-drug antibodies are a central readout in biotherapeutic development. An immune response raised against a molecule such as REGN5381 can alter its clearance, reduce exposure and blunt clinical response, so well-characterised ADA reagents are required as assay standards, positive controls and calibrators throughout preclinical and clinical immunogenicity testing.

Specifications

Product Name	Anti-REGN5381 Antibody (REGN9035)
SKU	GXDA00302
Target Drug	REGN5381
Antibody Type	Monoclonal
Clone ID	REGN9035
Host Species	Human
Reactivity	Human
Isotype	IgG
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Available Sizes	100 µg, 1 mg

Properties

Format	Liquid
Buffer	0.01M PBS, pH 7.4.
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Endotoxin Level	Please contact the lab for this information.
Guarantee	12 months from date of dispatch

Storage & Stability

Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.

Notes

For research use only. Not suitable for clinical or therapeutic use. Not for use in diagnostic or therapeutic procedures.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://www.assaygenie.com/content/Antibody%20Genie/GXDA/GXDA00302.pdf> | Version 1.0, 28 July 2026