

Research Grade Cendakimab

GXMA00012

Description

This Research Grade Cendakimab is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXMA00012
Contents:	Research Grade Cendakimab (50 µg, 100 µg) Bradford Reagent: 1 vial (2ml)
Target:	IMV surface fusion protein
Form:	Liquid
Storage Buffer:	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Product Usage Information:	ELISA: 1:2000-1:20000; WB: 1:500-1:2000

Antibody Data

Reactivity:	Monkeypox virus (MPXV)
Host species:	Mouse
Immunogen:	Recombinant MPXV A29L protein.
Clonality:	Monoclonal
Isotype:	IgG2b
Purification:	Protein A/G purified from cell culture supernatant.

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only.