

Research Grade Anrukinzumab

GXMA00017

Description

This Research Grade Anrukinzumab is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXMA00017
Contents:	Research Grade Anrukinzumab (50 µg, 100 µg) Bradford Reagent: 1 vial (2ml)
Target:	Thymidine kinase
Form:	Liquid
Storage Buffer:	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Product Usage Information:	ELISA, WB

Antibody Data

Reactivity:	African swine fever virus (ASFV)
Host species:	Mouse
Immunogen:	Recombinant ASFV K196R protein.
Clonality:	Monoclonal
Isotype:	IgG2a
Purification:	Protein A/G purified from cell culture supernatant.

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only.