

Anti-CSFV Envelope glycoprotein E1/gp33 Antibody (1A034)

GXMA00041

Description

This Anti-CSFV Envelope glycoprotein E1/gp33 Antibody (1A034) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXMA00041
Contents:	Anti-CSFV Envelope glycoprotein E1/gp33 Antibody (1A034) (50 µg, 100 µg) Bradford Reagent: 1 vial (2ml)
Target:	structural protein E1
Form:	Liquid
Storage Buffer:	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Product Usage Information:	ELISA, WB

Antibody Data

Reactivity:	Classical swine fever virus (CSFV)
Host species:	Mouse
Immunogen:	Synthesized peptide.
Clonality:	Monoclonal
Isotype:	IgG1
Purification:	Protein A/G purified from cell culture supernatant.

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only.