

Anti-Human Lambda light chain Monoclonal Antibody (1A181)

GXMA00165

Description

This Anti-Human Lambda light chain Monoclonal Antibody (1A181) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXMA00165
Contents:	Anti-Human Lambda light chain Monoclonal Antibody (1A181) (50 µg, 100 µg) Bradford Reagent: 1 vial (2ml)
Target:	Immunoglobulin lambda-1 light chain MCG
Form:	Liquid
Storage Buffer:	0.01M PBS pH 7.4
Product Usage Information:	ELISA, IHC, WB

Antibody Data

Reactivity:	Human
Host species:	Mouse
Immunogen:	Mammalian cells - derived recombinant human Lambda light chain protein.
Clonality:	Monoclonal

Isotype: IgG1

Purification: Protein A/G purified from cell culture supernatant.

Preparation & Storage

Stability and Storage: Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Store Bradford Reagent at Room Temperature for 1 Year.

Shipping: Gel Pack

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Note: For research use only.