

Anti-Human C6/Complement C6 Recombinant Nanobody (SAA1039)

GXRM01270

Description

This Anti-Human C6/Complement C6 Recombinant Nanobody (SAA1039) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXRM01270
Contents:	Anti-Human C6/Complement C6 Recombinant Nanobody (SAA1039) (50 µg, 100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	Complement component C6
Form:	Liquid
Storage Buffer:	0.01M PBS pH 7.4
Product Usage Information:	ELISA: 1:2000-1:20000

Antibody Data

Reactivity:	Human
Host species:	Alpaca
Clonality:	Monoclonal
Isotype:	VHH-8His-Cys-tag
Purification:	Purified by Nickel column.

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only.