

Human Rheumatoid Factor IgM Recombinant Monoclonal Antibody (2A2)

GXRM01840

Description

This Human Rheumatoid Factor IgM Recombinant Monoclonal Antibody (2A2) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXRM01840
Contents:	Human Rheumatoid Factor IgM Recombinant Monoclonal Antibody (2A2) (50 µg, 100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	Rheumatoid Factor
Form:	Liquid
Storage Buffer:	0.01M PBS pH 7.4
Product Usage Information:	Rheumatoid Factor (RF) Reference Standard

Antibody Data

Reactivity:	Human
Host species:	Human
Clonality:	Monoclonal
Isotype:	IgM
Purification:	Protein L purified from cell culture supernatant.

Preparation & Storage

Stability and Storage:	Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only.