

Anti-GIT1 Recombinant Monoclonal Antibody (N39B/8)

GXRM03485

Description

This Anti-GIT1 Recombinant Monoclonal Antibody (N39B/8) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXRM03485
Contents:	Anti-GIT1 Recombinant Monoclonal Antibody (N39B/8) (50 µg, 100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	CAT1
Form:	Liquid
Storage Buffer:	0.01M PBS pH 7.4
Product Usage Information:	WB: 1:500-1:2000

Antibody Data

Host species:	Mouse
Clonality:	Monoclonal
Isotype:	IgG2a, kappa
Purification:	Protein A/G purified from cell culture supernatant.

Preparation & Storage

Stability and Storage:	Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only.