

Anti-DENV-2 Envelope protein E/DIII C-C'loop Recombinant Monoclonal Antibody (E104#)

GXRM03671

Description

This Anti-DENV-2 Envelope protein E/DIII C-C'loop Recombinant Monoclonal Antibody (E104#) is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXRM03671
Contents:	Anti-DENV-2 Envelope protein E/DIII C-C'loop Recombinant Monoclonal Antibody (E104#) (50 µg, 100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Target:	Envelope protein E
Form:	Liquid
Storage Buffer:	0.01M PBS pH 7.4.
Product Usage Information:	ELISA: 1:2000-1:20000

Antibody Data

Reactivity:	Dengue virus 2 (DENV-2)
Host species:	Human
Clonality:	Monoclonal
Isotype:	IgG1, kappa

Purification: Protein A/G purified from cell culture supernatant.

Endotoxin Level: < 10 EU/mg

Preparation & Storage

Stability and Storage: Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Store Bradford Reagent at Room Temperature for 1 Year.

Shipping: Gel Pack

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Note: For research use only.