

Human IgM Recombinant Monoclonal Antibody

GXRM04216

Description

This Human IgM Recombinant Monoclonal Antibody is supplied for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	GXRM04216
Contents:	Human IgM Recombinant Monoclonal Antibody (50 µg, 100 µg, 1 mg) Bradford Reagent: 1 vial (2ml)
Form:	Liquid
Storage Buffer:	0.01M PBS pH 7.4

Antibody Data

Reactivity:	Human
Host species:	Human
Clonality:	Monoclonal
Isotype:	IgM
Purification:	Protein A/G purified from cell culture supernatant.

Preparation & Storage

Stability and Storage:	Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Gel Pack
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.
Note:	For research use only.