

Biotinylated Anti-equine IL-2 Polyclonal Antibody

KPAB0002

Description

This Biotinylated Anti-equine IL-2 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0002
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Goat
Reactivity: Equine
Immunogen: Recombinant Horse IL-2
Purification: Antigen-affinity Purification
Homology: Equus asinus (ass), Equus caballus (horse), Equus przewalskii (Przewalski's horse), Equus quagga (Plains zebra)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.