

## Anti-equine IL-4 Polyclonal Antibody

**KPAB0003**

### Description

---

This Anti-equine IL-4 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

### Product Information

---

|                       |  |
|-----------------------|--|
| <b>SKU:</b>           | KPAB0003                               |
| <b>Concentration:</b> | 1.0 mg/mL                              |
| <b>Conjugated:</b>    | Unconjugated                           |
| <b>Applications:</b>  | ELISA, ELISpot, FC, Neutralization, WB |

### Antibody Data

---

|                          |  |
|--------------------------|--|
| <b>Host:</b>             | Goat   |
| <b>Reactivity:</b>       | Equine   |
| <b>Immunogen:</b>        | Recombinant Horse IL-4   |
| <b>Molecular Weight:</b> | 12.8 kDa (calculated) for equine IL-4  |
| <b>Purification:</b>     | Antigen-affinity Purification  |
| <b>Formulation:</b>      | 1X PBS buffer containing 0.09% NaN <sub>3</sub>  |
| <b>Homology:</b>         | Equus asinus (ass), Equus caballus (horse), Equus przewalskii (Przewalski's horse), Equus zebra hartmannae (Hartmann's mountain zebra) |

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.