

## Biotinylated Anti-equine IL-4 Polyclonal Antibody

**KPAB0004**

### Description

---

This Biotinylated Anti-equine IL-4 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

### Product Information

---

**SKU:** KPAB0004

**Conjugated:** Biotin

**Applications:** ELISA, ELISpot

### Antibody Data

---

**Host:** Goat

**Reactivity:** Equine

**Immunogen:** Recombinant Horse IL-4

**Purification:** Antigen-affinity Purification

**Homology:** Equus asinus (ass), Equus caballus (horse), Equus przewalskii (Przewalski's horse), Equus zebra hartmannae (Hartmann's mountain zebra)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.