

## Anti-swine CCL2 Polyclonal Antibody

**KPAB0007**

### Description

---

This Anti-swine CCL2 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

### Product Information

---

<b>SKU:</b>	KPAB0007
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Swine
<b>Immunogen:</b>	Recombinant Porcine CCL2
<b>Molecular Weight:</b>	8.9 kDa (calculated) for swine CCL2
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X Tris buffer containing 0.09% NaN3
<b>Cross-Reactivity:</b>	Bovine CCL2: Strong; Canine CCL2: Strong; Equine CCL2: Strong; Guinea Pig CCL2: Weak; Mouse CCL2: Weak
<b>Homology:</b>	Sus scrofa (pig)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.