

## Anti-bovine CXCL9 Polyclonal Antibody

**KPAB0019**

### Description

---

This Anti-bovine CXCL9 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

### Product Information

---

<b>SKU:</b>	KPAB0019
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Bovine
<b>Immunogen:</b>	Recombinant Bovine CXCL9
<b>Molecular Weight:</b>	11.9 kDa (calculated) for bovine CXCL9
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN <sub>3</sub>
<b>Cross-Reactivity:</b>	Ovine CXCL9: Moderate
<b>Homology:</b>	Bison bison bison (American buffalo), Bos grunniens (domestic yak), Bos indicus (zebu), Bos mutus (wild yak), Bos taurus (cattle)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.