

Anti-swine IL-8 Polyclonal Antibody

KPAB0022

Description

This Anti-swine IL-8 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0022
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Swine
Immunogen:	Recombinant Porcine IL-8
Molecular Weight:	9.1 kDa (calculated) for swine IL-8
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN3
Cross-Reactivity:	Bovine IL-8: Strong; Canine IL-8: Moderate; Dolphin IL-8: Strong; Equine IL-8: Moderate; Feline IL-8: Strong; Ovine IL-8: Strong
Homology:	Phacochoerus africanus (Common warthog), Sus scrofa (pig)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.