

Anti-bovine IFN gamma Polyclonal Antibody

KPAB0025

Description

This Anti-bovine IFN gamma Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0025
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Bovine
Immunogen:	Recombinant Bovine IFN gamma
Molecular Weight:	16.9 kDa (calculated) for bovine IFN- η
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Canine IFN-H: Moderate; Caprine IFN-H: Strong; Dolphin IFN-H: Weak; Equine IFN-H: Weak; Feline IFN-H: Weak; Human IFN-H: Weak; Ovine IFN-H: Strong
Homology:	Bos frontalis (gayal), Bos taurus (cattle)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.