

Biotinylated Anti-canine IL-2 Polyclonal Antibody

KPAB0030

Description

This Biotinylated Anti-canine IL-2 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0030
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Canine
Immunogen: Recombinant Canine IL-2
Purification: Antigen-affinity Purification
Homology: Canis lupus dingo (dingo), Canis lupus familiaris (dog), Nyctereutes procyonoides (raccoon dog), Vulpes lagopus (Arctic fox), Vulpes vulpes (red fox)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.