

Anti-bovine IL-17A Polyclonal Antibody

KPAB0033

Description

This Anti-bovine IL-17A Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0033
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Bovine
Immunogen:	Recombinant Bovine IL-17A
Molecular Weight:	17.4 kDa (calculated) for bovine IL-17A
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.02% NaN ₃
Homology:	Bos indicus (zebu), Bos javanicus (banteng), Bos taurus (cattle)
Cross-Reactivity:	Canine IL-17A: Weak; Caprine IL-17A: Moderate; Cynomolgus IL-17A: Weak; Dolphin IL-17A: Strong; Feline IL-17A: Moderate; Ferret IL-17A: Weak; Guinea Pig IL-17A: Weak; Human IL-17A: Weak; Mouse IL-17A: Weak; Ovine IL-17A: Strong; Rabbit IL-17A: Weak; Rat IL-17A: Weak; Swine IL-17A: Strong

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.