

Anti-bovine TNF alpha Polyclonal Antibody

KPAB0034

Description

This Anti-bovine TNF alpha Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0034
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Bovine
Immunogen:	Recombinant Bovine TNF alpha
Molecular Weight:	17.4 kDa (calculated) for bovine TNF- α
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Ovine TNF-A: Moderate; Swine TNF-A: Weak
Homology:	Bison bison bison (American buffalo), Bos indicus (zebu), Bos javanicus (banteng), Bos mutus (wild yak), Bos taurus (cattle), Bubalus bubalis (water buffalo), Bubalus kerabau (carabao)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.